

LabKey User Conference 2019

Sample Lineage

Kevin Krouse
Susan Hert

October 4, 2019




The background of the slide is a solid blue color with a pattern of thin, light blue lines crisscrossing in various directions, creating a web-like or network effect.

What is Lineage?


Example: Sample Derivation



- Represents creation of daughter samples (e.g. aliquoting)
- Establishes hierarchy of data and samples along with metadata
- Linked to their parents via an experiment run
 - *(not shown on the graph)*




Example: Sample Derivation



Example: Sample Derivation




- Splitting
- Pooling
- Combined nodes for >5 siblings



Example: Assay + Transform (Summary)



Models assay process:
samples go in machine, file spit out, transformed file out




Example: Assay + Transform (Details)



Internal structure of
run captures the
transform step

Final outputs have
bold borders



Lineage is Experiment Runs



- Lineage backed by Experiment Run graph
- Experiment Run composition
 - Protocol Type - “Sample Derivation” or “MyAssay”
 - Inputs - Data or samples
 - Steps
 - Outputs - Data or samples

Experiment Run Use Cases



- Sample derivation
 - Splitting, pooling samples
 - Create parent/child lineage sample and data
- Assay Runs
- Pipeline Jobs sometimes
 - ETLs, Script pipelines
- XAR XML
 - Import/export Assay/SampleSet definitions, Runs


Creating Sample Lineage

Creating Sample Lineage



- Derive samples in LabKey Server UI
- Derive samples in Sample Manager UI
- Lineage columns for import/insert
 - Bulk import tsv
 - LABKEY.Query.insertRows
- LABKEY.Experiment.saveBatches/saveRuns
 - saveRuns - new in 19.3!

Deriving Samples: Sample Manager




Deriving Samples: Bulk Import



- Special columns for import
 - MaterialInputs/<SampleSetName>
 - DataInputs/<DataClassName>
- Comma separated list of names

Name	SampleDate	MaterialInputs/Samples	MaterialInputs/MixtureBatches	DataInputs/CellLine
SampleA	2019-09-10			
CellExample1	2019-09-10			CHO-K1
CellExample2	2019-09-10	CellExample1,SampleA		B-20190308-145

Deriving Samples: Bulk Import



B-20190308-145

Mixture Batches

[Overview](#) [Lineage](#)

Mixture Used	TBS
Expiration Date	2019-04-07
Consumed	false
Recipe Amount	8.0 L
Recipe Actual Amount	
Recipe Amount Unit	L

Ingredients Parents (3)

- [I-103](#)
- [I-58](#)
- [I-50](#)

Mixtures Parents (1)

- [TBS](#)

Samples Children (1)

- [CellExample2](#)

Deriving Samples: Query insertRows



```
LABKEY.Query.insertRows({  
  schemaName: "samples",  
  queryName: "Samples",  
  rows: [{  
    name: "Example1",  
    sampleDate: "2019-09-10",  
  }, {  
    name: "Example2",  
    sampleDate: "2019-09-10",  
    "MaterialInputs/Samples": "Example1,CellExample1"  
  }]  
});
```

Deriving: saveBatches/saveRuns API



```
LABKEY.Experiment.saveRuns({
  protocolName: LABKEY.Experiment.SAMPLE_DERIVATION_PROTOCOL,
  runs: [{
    name: "derive example",
    materialInputs: [{
      // or get existing sample by row id
      id: 105
    }, {
      // get existing sample by name
      name: "ES-1.205",
      sampleSet: {
        name: "ExpressionSystemSamples"
      }
    }
  ]},
  materialOutputs: [{
    // create a new sample
    name: "Example3",
    sampleSet: {
      name: "Samples"
    },
    properties: {
      SampleDate: "2019-09-10"
    },
  ]
}]);
```


Accessing Lineage

Accessing Lineage



- LabKey Server run graph UI
- LabKey Biologics / Sample Manager lineage UI
- LABKEY.Experiment.lineage() API
- Using Query to access lineage
- Using Lineage Query filter

Example: Lineage API



- Fetch lineage from a starting seed
- Indicate direction and depth
- Returns a graph of nodes

```
LABKEY.Experiment.lineage({  
  lsid: 'urn:lsid:labkey.com:Sample.31.Samples:KK01_t19119',  
  parents: true,  
  children: false,  
  depth: 3  
});
```

Example: Lineage API Response



```
{...}
└─ "urn:lsid:labkey.com:Sample.31.Samples:KK01_t19119": Start
  ├── children: []
  ├── cpasType: "urn:lsid:labkey.com:SampleSet.Folder-3235:Samples"
  ├── name: "KK01_t19119"
  └── parents: [ "urn:lsid:labkey.com:Run.Folder-3235:6031ac4f" ]
      ├── rowId: 1092270
      └── type: "Sample"
└─ "urn:lsid:labkey.com:Run.Folder-3235:6031ac4f":
  ├── children: [ "urn:lsid:labkey.com:Sample.31.Samples:KK01_t19119" ]
  ├── cpasType: "urn:lsid:labkey.org:Protocol:SampleDerivationProtocol"
  ├── name: "Derive 3 samples from KK01_t19114"
  └── parents: [ "urn:lsid:labkey.com:Sample.31.Samples:KK01_t19114" ]
      ├── rowId: 292262
      └── type: "Run"
└─ "urn:lsid:labkey.com:Sample.31.Samples:KK01_t19114":
  ├── children: [ "urn:lsid:labkey.com:Run.Folder-3235:6031ac4f" ]
  ├── cpasType: "urn:lsid:labkey.com:SampleSet.Folder-3235:Samples"
  ├── name: "KK01_t19114"
  └── parents: []
      ├── rowId: 1092265
      └── type: "Sample"
```

Lineage Query Columns



- Hidden Inputs and Outputs FK to lineage tables
 - Available on Run, Material, Data tables
 - Multi-valued foreign keys
- Performance implications
 - Lots of effort spent on lineage column perf
 - Lineage executed at render time
 - Sort/filter not possible

Example: Lineage Columns



Create New Update Status ▾ Assays Lineage ▾ ▾ ▾ + ▾ ▾ ▾ Derive Samples

Columns Filter Sort

Available Fields

- Outputs
 - All
 - Data
 - Materials
 - All
 - ☒ First
 - Samples
 - Runs

☒ Show Hidden Fields

Selected Fields

- Name
- Type
- Project
- Label
- Inputs/Materials/First
- Outputs/Materials/First

Delete Revert

View Grid Save

<input type="checkbox"/> ▾	Name	Type	Project	Label	Inputs/Materials/First	Outputs/Materials/First
<input type="checkbox"/>	KK01_t19116	Tissue	Kilo01	spleen	KK01_t19113	KK02_t19120, KK02_t19121
<input type="checkbox"/>	KK01_t19115	Tissue	Kilo01	spleen	KK01_t19113	KK02_t19120, KK02_t19121
<input type="checkbox"/>	KK01_t19113	Tissue	Kilo01	spleen		KK01_t19115, KK01_t19116

Lineage Query Filters



- LABKEY.Filter.Types.EXP_CHILD_OF
- Attach filter to an LSID column
- Filter value is LSID of a parent sample or data
 - `query.LSID~exp:childof=<parent-LSID>`

Example Application



demo time




- LabKey Documentation
 - [Experiment Terminology](#)
 - [Sample Derivation](#)
 - [Sample Status Demo](#) - ★ Premium Resource
- JS API docs
 - [LABKEY.Query.importData](#) - bulk import tsv
 - [LABKEY.Experiment.saveBatches](#)
 - [LABKEY.Experiment.lineage](#)

Extra Content

Example: XarTutorial Ex. 6 (Summary)




Run summary with inputs and outputs



Example: XarTutorial Ex. 6 (Details)




Lots of internal steps
with intermediate
samples and files



Example: Script Pipeline Run + Assay



RNASeqMatrixDemo



Example: Script Pipeline Run

